IntRaFast-Q SMA Tarama Kiti

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Spinal Müsküler Atrofi (SMA)

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Intelligent Ratio: A New Method for Carrier and Newborn Screening in Spinal Muscular Atrophy

Busranur Cavdarli,¹ Fatma Nihal Ozturk,² Sezen Guntekin Ergun,³ Mehmet Ali Ergun,⁴
Ozlem Dogan,⁵ and Emriye Ferda Percin⁴

Aim: Spinal muscular atrophy (SMA) is an inherited, autosomal recessive neuromuscular disease that causes high morbidity and mortality. The prevalence is 1–2/100,000, while the incidence is 1/6000–1/10,000 among live births. Due to the high carrier frequency (1/40–1/60) of SMA, screening can prevent new cases. The aim of the current study was to present the development of a new, quantitative, real-time, polymerase chain reaction (PCR)-based screening test that uses an intelligent ratio (IR) for analyses, as well as a comparison of the results with the gold standard.

Materials and Methods: Included in the study were 100 patients with various risk genotypes for survivor motor neuron 1 (*SMN1*) and *SMN2* genes whose genetics had been previously investigated using multiplex ligation probe amplification (MLPA). A combination of the 5' nuclease assay and allele-specific PCR was used to

Kaynaklar:

- ▶ 1. B. Cavdarli,1F. N. Ozturk, S. G. Ergun, M. A. Ergun, O. Dogan and E. F. Percin. "Intelligent Ratio:A New Method for Carrier and Newborn Screening in Spinal Muscular Atrophy". Genetic Testing And Molecular Biomarkers, Volume 24, Number 9, 2020.
- ▶ 2. Shin S1, Park SS, Hwang YS, Lee KW, Chung SG, Lee YJ, Park MH. Deletion of SMN and NAIP genes in Korean patients with spinal muscular atrophy. J Korean Med Sci. 2000;15:93-8.
- ▶ 3. Verhaart IEC, Robertson A, Wilson IJ, Aartsma-Rus A, Cameron S, Jones CC, Cook SF, Lochmüller H. Prevalence, incidence and carrier frequency of 5q-linked spinal muscular atrophy - a literature review. Orphanet J Rare Dis. 2017;12:124.
- ▶ 4. Wirth B. An update of the mutation spectrum of the survival motor neuron gene (SMN1) in autosomal recessive spinal muscular atrophy. Hum Mutat. 2000;15:228–37.
- ▶ 5. Boardman FK, Young PJ, Warren O, Griffiths FE. The role of experiential knowledge within attitudes towards genetic carrier screening: A comparison of people with and without experience of spinal muscular atrophy. Health Expect. 2017 Jul 13. doi: 10.1111/hex.12602.
- ▶ 6. Shuji Ogino, Debra G. B. Leonard, Hanna Rennert, and Robert B. Wilson. Spinal Muscular Atrophy Genetic Testing Experience at an Academic Medical Center. J Mol Diagn. 2002 Feb; 4(1): 53–58.
- ▶ 7. Han C. Phan, Jennifer L. Taylor, Harry Hannon and Rodney Howel. Newborn screening for spinal muscular atrophy: Anticipating an imminent need. Seminars in Perinatology. April 2015 Volume 39, Issue 3, Pages 217–229.
- ▶ 8. Thomas W Prior, Narasimhan Nagan, Elaine A Sugarman, Sat Dev Batish and Corey Braastad. Technical standards and guidelines for spinal muscular atrophy testing. Genetics in Medicine volume 13, pages 686–694 (2011)
- ▶ 9. Jin He, Qi-Jie Zhang, Qi-Fang Lin, Ya-Fang Chen, Xiao-Zhen Lin, Min-Ting Lin, Shen-Xing Murong, Ning Wang and Wan-Jin Chen. Molecular analysis of SMN1, SMN2, NAIP, GTF2H2, and H4F5 genes in 157 Chinese patients with spinal muscular atrophy. Gene Volume 518, Issue 2, 15 April 2013, Pages 325-329
- ▶ 10. Amara A1, Adala L, Ben Charfeddine I, Mamaï O, Mili A, Lazreg TB, H'mida D, Amri F, Salem N, Boughammura L, Saad A and Gribaa M.. Correlation of SMN2, NAIP, p44, H4F5 and Occludin genes copy number with spinal muscular atrophy phenotype in Tunisian patients.Eur J Paediatr Neurol. 2012 Mar;16(2):167-74
- ▶ 11. American College of Obstetricians and Gynecologists' Committee on Genetics in collaboration with committee members Britton Rink, Stephanie Romero, Joseph R Biggio Jr, Devereux N. Saller Jr. and Rose Giardine.Committee Opinion. Number 691. March 2017.

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© Mary Ann Liebert, Inc.